日本加速器学会年会、幕張、2016年8月8-10日

MeVフェムト秒電子線パルスを用いた 超高速電子顕微鏡の開発

楊 金峰

大阪大学·産業科学研究所







Motivation

✓An observation in "real space(実空間)" and "real time(実時間)" has long been in goal for the scientists.

 Many ultrafast phenomena or dynamics in materials begin in femtosecond time region over nanometer spatial dimension, e.g.

- Ultrafast chemical reactions (e.g. excitation, relaxation, making and breaking of bonds, ...)
- Ultrafast dynamics in biology (e.g. charge/energy transfers, growth of damage, ...)
- Structural changes in solids (e.g. phase transformations, melting and resolidification, ...)
- ✓ Single-shot measurement is necessary, because many structural phenomena or dynamics are not reversible (不可逆的な構造変化)!

MeVフェムト秒電子線パルスを用いた超高速電子顕微鏡の開発

フォトカソードRF電子銃



Concept of RF gun based UEM



Femtosecond photocathode electron gun



Electron energy :	1 ~ 3 MeV
Bunch length :	100 fs
Emittance :	0.1 mm-mrad
Energy spread :	10 ⁻⁴ ~ 10 ⁻⁵
Charge:	10 ⁷ ~10 ⁸ e⁻'s/pulse

Time resolution:100 fsSpatial resolution:1 nm



RF電子銃を用いた 超高速電子顕微鏡実証機の写真



$0.1 \mu m$ - emittance electron beam generation

Measurement results of emittance

(1) <u>reduce the thermal emittance</u> by focusing laser on the cathode



for emittance measurement

Dependence of image contrast on emittance



フェムト秒電子線パルスによる電子回折の測定結果





<u>Sample:</u> Single-crystal Au Thickness: ~10nm

Electron beam: Energy: 3.1 MeV Pulse length: 100 fs e- charge: ~1pC/pulse (10⁷e-'s/pulse)

Results:

- High-quality MeV electron diffraction patterns can be observed with MeV-highenergy e- pulses.
- Single-shot measurement is available.

MeV TEM imaging of MoO₃ crystal

<u>Electron beam</u>: 3.1 MeV, 100 fs, ~1pC/pulse Diameter of MoO₃ crystal: ~2μm





Dependence on electron pulse number

<u>Electron beam</u>: 3.1 MeV, 100 fs, ~1pC/pulse <u>Sample</u>: MoO₃ crystal (diameter: ~2 μ m)



➤Contrast TEM images can be observed with 500 pulses.

シングルショットとパルス積算で観測したTEM像の比較

<u>Electron beam</u>: 3.1 MeV, 100 fs, ~1pC/pulse <u>Sample</u>: single crystal gold (thickness: 10nm)



- In the averaging measurement, good-contrast TEM images can be observed for Au sample.
- >At low-magnification observation, the single-shot imaging is available.

まとめ

- ✓ 世界先駆けて相対論的フェムト秒電子線パルスを用いた超高速電子顕 微鏡実証機を開発し、原理実証に成功した。
- ✓ 電子回折の測定モードでは、シングルショットの観測でも十分明瞭な電子 回折パターンを得ることができる。
- ✓ 電子顕微鏡像の測定モードでは、3MeVのフェムト秒電子線パルスを用いて、微粒子や微結晶等のTEM像の観測に成功し、低倍率ではシングルショットの観測も可能である。

<u>今後の課題:</u>

- ✓ 電子ビーム高輝度化 → エミッタンスの低減(<0.1µm)、エネルギー分散の低減(~10⁻⁵)
- ✓ 安定度(~10⁻⁵)と検出効率の向上
 - ⇒ 空間分解能の向上、フェムト秒超高速電子顕微鏡の実現!

謝 辞

<u>共同研究者</u>

 (阪大産研)淺川稜、菅晃一、近藤孝文、神戸正雄、吉田陽一 (阪大産研)吉田秀人、竹田精治 (阪大工)谷村克己、保田英洋 (北大)成瀬延康 (ローザンヌ工科大)室岡義栄 (KEK)高富俊和、照沼信浩、浦川順治

<u>産研量子ビーム科学研究施設のご協力</u>



基盤研究(A)(22246127、H22~24)、基盤研究(A)(26246026、H26~28)、 萌芽研究(16K13687、H28~29)